

# Sapphire Fine-Light Flexi. The most amazing product

The Sapphire Fine-Light Flexi is a safe and revolutionary skin rejuvenation system combining the latest Red LED Light Therapy technology with an active Hyaluronic Acid/Collagen Booster mask.

The Sapphire Fine-Light Flexi mask is easily portable, and in just minutes combined with our Fine-Light Flexi Collagen Booster will assist your skin to display a more youthful look. You will notice the difference in no time at all.

The Sapphire Fine-Light Flexi is light, comfortable and easily worn.

### LED Light Therapy - how it works?

Red LED light (635nm) is a powerful stimulant to the dermis and is an effective way to energise the cells of the skin, promoting the production of collagen and elastin. Lack of collagen and elastin in the skin's cells can result in fine lines and wrinkles. LED Light Therapy has a direct effect on the regeneration of the skin, penetrating the top layers of the skin, energising the skin, producing collagen and elastin, which helps to plump the skin and reduce fine lines and wrinkles. This is known as Photo Rejuvenation.

## **Light Therapy Works!**

Light Therapy has the desired effect of rejuvenating the skin without the need for surgery, injections or painful procedures. Referred to as the 'No Knife Face Lift' this technology is fast becoming the preferred solution for those seeking to strip away the years. It is the safest, most affordable way to achieve radiant, younger looking skin.

So, why can't you get the same effect from sunlight? Well, you can, but sunlight also contains harmful Ultra Violet (UV) rays which damage the skin. Fine-Light Flexi bathes the skin in a pure light source that generates much needed protein and in turn, allows the cells of the dermis to regenerate, plumping the skin and increasing levels of elasticity.

#### Does LED Light Therapy really work?

LED Light Therapy technology is tried and tested. Used extensively in professional salons, this technology is now available for home use. Tested and evaluated by many medical practitioners, the technology was developed by NASA in the United States and continues to be the focus of reports and white papers on skin rejuvenation.

### What is the active Hyaluronic Acid/Collagen Booster mask?

The Sapphire Fine-Light Flexi Collagen Mask is placed directly on the face under the Sapphire Fine-Light Flexi Red Light Mask. Active ingredients, such as high dosed collagen and hyaluronic acid promote regeneration and give the skin new radiance.

#### Its so easy to use

The Fine-Light Flexi combination therapy is so easy to use. Simply place the unique collagen booster and red light mask directly over the face, turn it on and leave for 20 minutes while your client relaxes.

Is anti-aging LED Light Therapy right for you?

Of course it is! Control the aging with the Sapphire Fine-Light Flexi.



Laser SOS Aesthetics Ltd. Unit 3 Burrel Road, St Ives, Cambs PE27 3LE United Kingdom

Tel: +44 (0)1480 409962 Fax: +44 (0)1480 498860

e-Mail: sales@lasersosaesthetics.com

www.lasersosaesthetics.com



A World Class Service from a World Class Company